

SARS-CoV-2 Environmental Sampling and qPCR Analysis

Jason Dobranic, Ph.D., CEICC, CIEC
Vice President of Microbiology
& Life Sciences at EMSL Analytical, Inc.

EMSL

SM





ASBESTOS DIVISION

^{82}Pb

LEAD DIVISION



MICROBIOLOGY DIVISION



INDUSTRIAL HYGIENE DIVISION



ENVIRONMENTAL CHEMISTRY DIVISION



MATERIALS SCIENCE DIVISION




RADIOCHEMISTRY DIVISION



FOOD AND BEVERAGE DIVISION

EMSL Analytical, Inc. is one of the leading testing laboratories with 46 locations throughout the United States and Canada.

United States and Canada Locations





Transmission and Surface Contamination

- Main infection route is between people who are in close contact (within about 6 feet).
- Through respiratory droplets produced when an infected person coughs or sneezes.
- Currently, transmission from surfaces is not thought to be the main way the virus spreads, however the CDC's April 13th update remarked:
"COVID-19 is a new disease and we are still learning about how it spreads and the severity of illness it causes."
- Possible by touching a surface or object that has the virus on it, and then touching the mouth, nose, or eyes.



Environmental Surfaces

- May remain infective on surfaces for hours to days
- Dependent on the surface material type (fabric, tile, steel, paper, etc.)
- Thorough cleaning and disinfection of frequently touched surfaces are recommended by the CDC and believed to be essential in preventing the spread of infection





Stability of SARS-CoV-2 in Different Environmental Conditions

- Virus is highly stable at 4°C (14 days), but sensitive to heat
- Temperature increased to 70°C, the time for virus inactivation was reduced to 5 mins
- No infectious virus could be recovered from printing and tissue papers after a 3-hour incubation
- Day 2 treated wood and cloth
- Day 4 glass and banknote
- Day 7 stainless steel and plastic
- Surgical mask on day 7





Potential High Risk Areas to Sample

- Telephones, Keyboards
- Chairs, Tables, Light Switch
- Hand Soap Dispenser
- Garbage Bin
- Office Equipment, Panels
- Doorknobs
- Railings, banisters
- Faucet Handles
- Toilet Handles
- Air Returns, Ventilation Exits





Disinfectants Studied

Disinfectant (Working concentration)	Virus titre (Log TCID ₅₀ /mL)		
	5 mins	15 mins	30 mins
Household bleach (1:49)	U	U	U
Household bleach (1:99)	U	U	U
Hand soap solution (1:49)	3.6 [#]	U	U
Ethanol (70%)	U	U	U
Povidone-iodine (7.5%)	U	U	U
Chloroxylonol (0.05%)	U	U	U
Chlorhexidine (0.05%)	U	U	U
Benzalkonium chloride (0.1%)	U	U	U



ATP Testing


- Adenosine Triphosphate – the energy molecule used by living cells
- Measured by luminometer (RLUs) onsite
- Commonly used in food processing
- Indicator of cleanliness (all cellular ATP)
 - food residues, blood spillages, fecal matter, dust, shed skin squama
- Not produced or stored by viruses
- Not a direct measure of viral presence
- Follow manufacture recommendation on sampling





Microbial Screening

- Swab surfaces for bacteria/fungi as indicators of cleaning efficacy
- Assumes an EPA-registered virucide was included in the cleaning protocol and used according to labeling
- If results are free of bacteria and fungi then those surface have been properly cleaned




<https://www.worldbioproducts.com/purblue.html>



Microbial Screening

- Wait until disinfectant has dried
- Must use swab with neutralizing buffer
- Fungal plate counts, 5 day TAT
- Bacterial plate counts, 2 day TAT
- Ship swabs on ice packs



<https://www.worldbioproducts.com/purblue.html>



SARS-CoV-2 PCR Testing

- Presence or absence of the viral genetic material on environmental surfaces
- RT-qPCR targets the SARS-CoV-2-nucleocapsid N gene
- Two indicators: 2019-nCoV_N1 & 2019-nCoV_N2
- Same method as specified by the CDC and authorized by the FDA for clinical COVID-19 diagnostics
- Human RNase P gene is not analyzed
- Not for processing patient SARS-CoV-2 specimens



SARS-CoV-2 PCR Testing cont.

- Environmental samples frequently contain qPCR inhibiting substances found in dirt, soil, and dust
- Additional nucleic acid (NA) extraction method (clean-up) for environmental samples
- 'Detected' for either indicator can be interpreted with confidence that the area swabbed contained SARS-CoV-2 coronavirus particles that could be deemed as potentially infectious
- Test includes five controls: positive N1, positive N2, negative N1, negative N2, and qPCR inhibition control





SARS-CoV-2 PCR Results

- Test is currently P/A results but may be quantitative in future
- Limit of detection (LOD) 1 RNA copy per microliter of transport medium
- The test is specific only to SARS-CoV-2
- There is no known cross-reaction to any other viruses
- Rapid TATs available: same day to 3 day typically offered by commercial labs



Swab Sampling Protocol

- Wear appropriate PPE during sampling
- Swab surface areas of interest (25 cm² is commonly used)
- Moderate pressure while moving in at least two different directions and rotating the swab
- Avoid letting the swab dry completely





Swab Sampling Protocol (cont.)

- Place back in labeled transport container
- Wipe with disinfectant
- Place all swabs in Ziplock bag, wipe bag
- Double bag and wipe outside of second bag
- Place COC in separate Ziplock bag





Air Sampling

- Vira-Pore - 37mm 1 μm pore-size PTFE-laminated PTFE Zefon[®] "ZePore"[™] filters
- Tested with 10 L/min flow rate for 10 mins
- Other samplers
 - BioSampler, Button Sampler (gelatin, short time)



SKC

Specialists in Air Sampling Since 1962




**ENVIRONMENTAL
EXPRESS[®]**



VIRA-PORE Viral Sampling Cassette

COVID-19 in Wastewater





Wastewater Sampling

- Collect wastewater as it enters treatment plant (untreated influent)
- Wear gloves and do not touch the bottle opening
- Sterile 120mL bottle (Only 105ml required for analysis)
- Wipe outside of the bottles, then double bag
- Ship same day to the lab in cooler with ice packs
- Samples may be stored up to 72 hours at 2-8°C



Wastewater Samples

- IDEXX Water SARS-CoV-2 RT-PCR Test Kit (P/A)
 - Detects the same virus regions targeted during the analysis of environmental swabs or air samples
 - Sterile 120mL Bottle (Only 105ml required for analysis)
- New tests being developed for quantitative results





Sampling Strategy

- Detect the presence of COVID-19 within a sewershed
- Monitor trends in COVID-19 within a sewershed
- Targeted wastewater surveillance entails sampling wastewater from upstream in the wastewater network (e.g., lift stations, interceptors, manholes)
- SARS-CoV-2 RNA concentrations in wastewater cannot currently be used to determine the total number of infected persons in a community or the percent of the population that is infected.



Sampling Strategy

- A minimum of three samples is required to detect wastewater trends over time.
- Grab samples can be collected rapidly but are less representative
- Composite samples are collected by pooling multiple grab samples at a specified frequency over a set time period – typically 24 hours for wastewater surveillance.

From CDC - <https://www.cdc.gov/coronavirus/2019-ncov/cases-updates/wastewater-surveillance/developing-a-wastewater-surveillance-sampling-strategy.html>



Shipping Samples

- Pack well so samples are not damaged during shipping
- Samples must be shipped with an ice-pack
- Recommend shipping samples on the same day as collected
- Samples may be refrigerated up to 72 hours per WHO guidelines





Interpreting Results

- Results are reported in either 'Detected' or 'Non-Detected' for two indicators of the SARS-CoV-2
- Indicator one is 2019-nCoV_N1 and indicator two is 2019-nCoV_N2.
- Both indicators detect the nucleocapsid protein 'N' gene of the SARS-CoV-2 coronavirus.
- Either indicator can be interpreted that the area sampled did contain SARS-CoV-2 coronavirus particles that could be deemed as potentially infectious.



EMSL ANALYTICAL, INC.
TESTING LABS • PRODUCTS • TRAINING

Jason Dobranic, Ph.D.

800-220-3675

jdobranic@emsl.com

www.emsl.com

